


# SpectraMax Multi-Mode Microplate Readers

YOUR APPLICATIONS, YOUR MODES, YOUR CHOICE



- UPGRADEABLE PLATFORM FOR CHANGING LAB NEEDS
- THREE-MODE CUVETTE PORT FOR ASSAY DEVELOPMENT
- DUAL MONOCHROMATOR TUNABILITY
- AUTOMATED ABSORBANCE PATHLENGTH CORRECTION
- ENDPOINT, KINETIC, SPECTRAL AND WELL-SCANNING READ TYPES
- COMPREHENSIVE DATA ANALYSIS WITH SOFTMAX PRO SOFTWARE
- VALIDATION & COMPLIANCE
- ROBOTICS COMPATIBILITY

The SpectraMax® M3, M4, M5 and M5<sup>c</sup> Multi-Mode Microplate Readers are a modular, upgradeable dual-monochromator microplate reader platform offering a wide range of high performance multi-mode capabilities ideal for life science research and drug discovery screening. Choose from a three- (M3), four- (M4), or five- (M5/M5<sup>c</sup>) mode reader customized to your specific applications or budgetary needs, while optional capabilities allow you to upgrade with other detection modes at a later time. All configurations offer a triple-mode cuvette port, accurate temperature control, microplate shaking and comprehensive data management using our SoftMax® Pro Software. Detection modes include:

- UV-Visible Absorbance (Abs)
- Fluorescence Intensity (FI)
- Luminescence (Lum)
- Time-Resolved Fluorescence (TRF)
- Fluorescence Polarization (FP)

The SpectraMax M5<sup>c</sup> Reader offers the additional benefit of being certified for Cisbio Bioassays' HTRF® Assays.


## DUAL MONOCHROMATORS FOR ASSAY FLEXIBILITY

With SpectraMax Multi-Mode Readers, there is no need to utilize expensive filters to optimize detection levels and background. The optical systems use two scanning monochromators so the user can determine optimal excitation and emission settings, resulting in assay performance similar to that of dedicated single-mode readers.

## PATENTED PATHLENGTH CORRECTION FOR BETTER ABSORBANCE ACCURACY

Only Molecular Devices microplate plate readers offer the capability to measure the depth (optical pathlength) of samples with no temperature dependency using the patented PathCheck® Sensor technology. With SoftMax Pro Software, the PathCheck Sensor automatically normalizes the well absorbance. This eliminates the need for standard curves, and, for compounds with known absorptive properties, enables users to calculate concentrations directly from absorbance.

## Superior Optics for Optimal Assay Performance




### UNIQUE OPTICAL CHARACTERISTICS

1. The reference diodes enable elimination of measurement noise due to slight fluctuations in excitation light intensity.
2. The angled emission beam improves signal-to-noise, especially in narrow Stokes shift fluorophores, by reducing stray light.
3. Elliptical mirrors are used instead of lenses for maximum transmission with minimal wavelength distortion.
4. Top-quality UV-grade fibers give the highest light transmission down to even the lowest wavelengths.

### ASSAY COLLABORATION FOR EASE OF SETUP

Molecular Devices has collaborated with various assay partners to optimize and validate homogeneous and heterogeneous biochemical- or cell-based assay performance on the SpectraMax platform. To support these assays, we provide application notes as well as ready-to-run protocols in our SoftMax Pro Software. Some of our featured partner assays include HTRF assays from Cisbio Bioassays and Lanthascreen® TR-FRET assays from Invitrogen.




## Which SpectraMax Microplate Reader Do You Need?

	SpectraMax M3 Reader	SpectraMax M4 Reader	SpectraMax M5 Reader	SpectraMax M5 <sup>e</sup> Reader
<b>Detection Modes</b>				
Absorbance	✓	✓	✓	✓
Fluorescence Intensity	✓	✓	✓	✓
Luminescence	✓	✓	✓	✓
Time-Resolved Fluorescence		✓	✓	✓
Fluorescence Polarization			✓	✓
HTRF				✓
Upgrade Options	TR,F,HTRF,FP	HTRF, FP	HTRF	N/A
<b>Plate Formats</b>				
6-, 12-, 24-, 48-, 96-, 384-Well Microplates	✓	✓	✓	✓
<b>Certification &amp; Validation</b>				
IMAP Validation		✓ (TR-FRET only)	✓ (TR-FRET/FP)	✓ (TR-FRET/FP)
HTRF Certification				✓
LanthaScreen Certification			✓	✓
<b>Key Applications</b>				
ELISAs	✓	✓	✓	✓
DNA/RNA Quantitation	✓	✓	✓	✓
Protein Assays	✓	✓	✓	✓
Enzyme Kinetics	✓	✓	✓	✓
Protease Assays	✓	✓	✓	✓
Fluorescent Proteins and FRET	✓	✓	✓	✓
Cell Viability and Cytotoxicity Assays	✓	✓	✓	✓
Reporter Gene Assays	✓	✓	✓	✓
Cell Migration Assays	✓	✓	✓	✓
QBT Fatty Acid Uptake Assay	✓	✓	✓	✓
Neurotransmitter Transporter Uptake Assay	✓	✓	✓	✓
ADME-Tox	✓	✓	✓	✓
Membrane Permeability	✓	✓	✓	✓

## Software for Data Acquisition and Analysis

### Powerful, Easy-to-Use Software



SoftMax Pro Software features easy-to-set-up protocols for data acquisition, customizable spreadsheet functionality for analysis, and powerful graphing tools for data presentation.

### COMPREHENSIVE DATA ANALYSIS AND GxP SOLUTIONS

SoftMax Pro Software provides data acquisition, analysis and management capabilities, allowing cross-plate analysis and custom calculations. For users operating in a FDA 21 CFR Part 11 compliant environment, SoftMax Pro GxP Software is available, allowing user permissions, audit trails, e-signature and reporting tools.

## Tools for Validation and Compliance

### Validation Test Plates for Abs, FL, Lum Optical Performance



SpectraTest ABS1 Absorbance Plates, FL1 Fluorescence Plates, and LM1 Luminescence Plates are used to validate optical performance of SpectraMax M3, M4, and M5/M5e Microplate Readers.

### VALIDATION AND COMPLIANCE OF OPTICAL CHARACTERISTICS

SpectraMax Readers have the most complete level of product validation and compliance. Molecular Devices provides the complete solution covering the instrument and software:

- SpectraTest® ABS1, FL1, and LM1 Validation Plates for hardware validation of absorbance, fluorescence, and luminescence modes
- IQ/OQ/PQ for all microplate readers
- SoftMax Pro Software Validation Package
- Software tools for FDA 21 CFR Part 11 compliance

## Automation Solutions

### StakMax Microplate Handling System



Integrate any SpectraMax Multi-Mode Reader with the StakMax Microplate Handling System from Molecular Devices. The system provides automation for up to 50 microplates for easy walkaway automation. System setup and calibration are controlled from within SoftMax Pro Software.

### Compatible Automation Solutions for SpectraMax Readers



SoftMax Pro 5 Software has been integrated by many leading robotics and LIMS providers, enabling both data analysis and instrument control in automated environments.

### ROBOTICS COMPATIBILITY FOR INCREASED THROUGHPUT

SpectraMax Multi-Mode Microplate Readers can be easily integrated with our optional StakMax<sup>®</sup> Microplate Handling System for walk-away processing. Operated from within SoftMax Pro Software, the StakMax System can hold up to 50 plates and facilitates barcode reading.

For more advanced automation needs Molecular Devices interacts with all of the major lab automation providers, and is one of their leading choices.

### ORDERING INFORMATION

Acquiring a SpectraMax Multi-Mode Microplate Reader is extremely easy:

1. Decide what modes you need and choose the specific configuration option that suits you best. All systems include SoftMax Pro Software for Windows<sup>®</sup> and Macintosh<sup>®</sup> Operating Systems.
2. Choose the additional options you want:
  - Software validation tools
  - SpectraTest ABS1, FL1, LM1 Validation Plates
  - SoftMax Pro GxP Software
  - StakMax Microplate Handling System
3. Contact your Molecular Devices sales representative to discuss the details.



## Technical Specifications

### General Specifications

Dimensions (in.): 8.6 (H) x 22.8 (W) x 15.3 (D)  
 Dimensions (cm): 22 (H) x 58 (W) x 39 (D)  
 Weight: 36 lbs. (16.4 kg)  
 Power consumption: < 420 watts  
 Power source: 100–240 VAC, 3.5 A, 50/60 Hz  
 Robotic-compatible: Yes

### General Photometric Performance

Plate formats: 6, 12, 24, 48, 96, 384 wells  
 Light source: Xenon Flash Lamp  
 (1 joule/flash)  
 Detectors: 2 photomultiplier tubes (PMT)  
 Shaker time: 0 to 999 seconds  
 Temp. control: 2°C above ambient to 60°C  
 Temp. uniformity: < 1°C at 37°C set point  
 Temp. accuracy: ±1°C at 37°C set point  
 Endpoint reading: All modes  
 Kinetic reading: All modes  
 Spectral scanning: All modes  
 Well scanning: Abs, FI, TRF, Lum

### Typical Read Times (minutes:seconds)\*

	96 wells	384 wells
Absorbance	0:18	0:49
Fluorescence Intensity	0:17	0:48
Fluorescence Polarization	0:42	2:03
Time-Resolved Fluorescence	0:17	0:48
Luminescence	2:00	7:00

\*With 3 flashes/well in absorbance and fluorescence modes, and 1 sec./well integration in luminescence.

### Absorbance Photometric Performance

Reading capabilities: Cuvette or microplate  
 Wavelength range: 200–1000 nm  
 Wavelength selection: Monochromator, tunable in 1.0 nm increments  
 Wavelength bandwidth: ≤ 4.0 nm  
 Wavelength accuracy: ±2.0 nm  
 Wavelength repeatability: ±0.2 nm  
 Photometric range: 0–4.0 OD  
 Photometric resolution: 0.001 OD  
 Photometric accuracy (microplate): < ±0.006 OD ±1.0%, 0–2 OD  
 Photometric accuracy (cuvette): < ±0.005 OD ±1.0%, 0–2 OD  
 Photometric precision: < ±0.003 OD ±1.0%, 0–2 OD  
 Stray light: < 0.05% @ 230 nm

### Fluorescence Intensity Performance

Reading capabilities: Cuvette or top or bottom of a microplate  
 Wavelength range: 250–850 nm  
 Wavelength selection: Monochromators, tunable in 1.0 nm increments  
 Bandwidth (EX, EM): 9 nm, 15 nm  
 Sensitivity: < 5 pM fluorescein in 96 wells or cuvette, < 20 pM in 384 wells

### Luminescence Performance

Reading capabilities: Cuvette or top or bottom of a microplate  
 Wavelength selection: Choice of simultaneous detection of all wavelengths or selection via monochromator, tunable in 1.0 nm increments  
 Wavelength range: 250–850 nm  
 Sensitivity: < 2 fg/well lower detection limit for firefly luciferase in 96- and 384-well top read  
 Dynamic range: > 6 decades  
 Cross-talk: < 0.3% in white 96- and 384-well microplates

### Time-Resolved Fluorescence Performance (M4, M5, M5<sup>e</sup> only)

Reading capabilities: Top or bottom of a microplate  
 Wavelength range: 250–850 nm  
 Wavelength selection: Monochromators, tunable in 1.0 nm increments  
 Bandwidth (EX, EM): 9 nm, 15 nm  
 Precision data collection: 1–100 flashes, delay of 0–600 μsec. before read, integration time selectable between 50–1500 μsec.  
 Sensitivity: 100 fM europium in 96 or 384 wells with top-read  
 SpectraMax M5<sup>e</sup> reader only: Certified to Cisbio Bioassays' HTRF assays performance specifications

### Fluorescence Polarization Performance (M5/M5<sup>e</sup> only)

Wavelength range: 400–750 nm  
 Wavelength selection: Monochromators, tunable in 1.0 nm increments  
 Bandwidth (EX, EM): 9 nm, 15 nm  
 Precision: < 5 mP standard deviation at 1 nM fluorescein in 96 and 384 wells

### Patents

The PathCheck Sensor is covered under U.S. Patents 5,959,738, 6,188,476, 6,320,662, 6,339,472, 6,404,501, 6,496,260 and 6,995,844. SpectraMax M3, M4, M5, and M5<sup>e</sup> readers are also covered under U.S. Patents 6,097,025, 6,232,608, 6,236,456, 6,313,471 and 6,316,774, 6,693,709, and 6,825,921.

### SALES OFFICES

#### United States & Canada

Molecular Devices  
 Tel. +1-800-635-5577  
 Fax +1-408-747-3601

#### Brazil

Molecular Devices Brazil  
 Tel. +55-11-3616-6607  
 Fax +55-11-3616-6607

#### China

Molecular Devices Beijing  
 Tel. +86-10-6410-8669  
 Fax +86-10-6410-8601

#### Molecular Devices Shanghai

Tel. +86-21-6887-8820  
 Fax +86-21-6887-8890

#### Germany

Molecular Devices GmbH  
 Tel. +49-89/96-05-88-0  
 Fax +49-89/96-02-34-5

#### Japan

Molecular Devices Japan, Osaka  
 Tel. +81-6-6399-8211  
 Fax +81-6-6399-8212

#### Molecular Devices Japan, Tokyo

Tel. +81-3-5282-5261  
 Fax +81-3-5282-5262

#### South Korea

Molecular Devices Korea, LLC  
 Tel. +82-2-3471-9531  
 Fax +82-2-3471-9532

#### United Kingdom

Molecular Devices (GB) Ltd.  
 Tel. +44-118-944-8000  
 Fax +44-118-944-8001

[www.moleculardevices.com](http://www.moleculardevices.com)

FOR RESEARCH USE ONLY. NOT FOR USE IN  
 DIAGNOSTIC PROCEDURES.

The trademarks used herein are the property of  
 Molecular Devices, Inc. or their respective owners.

Specifications subject to change without notice.