

Gold Standard TaqMan® Assay Technology Enables Highly Accurate Quantitation for Copy Number Variation

TaqMan® Copy Number Assays are the gold standard for accurate target quantitation, making them ideal for validating and screening specific copy number variations. These assays combine TaqMan® Assay probe 5' nuclease chemistry with Applied Biosystems® real-time PCR instruments to form a method for obtaining specific, reproducible, and easy-to-interpret results. TaqMan® Copy Number Assays are run together with a TaqMan® Copy Number Reference Assay in a duplex real-time polymerase chain reaction (PCR). Results are analyzed by the relative quantitation method using CopyCaller™ Software. On average, setup to primary analysis takes 3 to 4 hours, enabling users to obtain results with speed, confidence, and reproducibility.

Features and Benefits

1. Performance

Accuracy

Assays are developed with a high degree of bioinformatics to overcome homology issues and help ensure accurate detection of a wide range of specific variations. Gain confidence in studying the right copy number variation to generate correct results. CopyCaller™ Software provides the calculated copy number and predicted copy number, along with confidence value and z-score quality metrics.

Specificity

TaqMan® Copy Number Assays use TaqMan® Assay chemistry and Applied Biosystems® Real-Time PCR Systems to obtain specific and reproducible copy number results.

Minor Groove Binder (MGB) moiety (5' FAM® reporter dye/3' MGB/nonfluorescent quencher (NFQ)):

- Stabilizes the hybridized probe, resulting in stronger probe binding
- Allows the probe to be shorter than traditional dual-labeled probes, enabling you to effectively target difficult sequences and increase specificity

2. Design and Coverage

TaqMan® Copy Number Assays are designed using Applied Biosystems® assay design algorithms, which are optimized to produce high-performing copy number assays. TaqMan® Copy Number Assays are based on the NCBI Reference Genome Build 37.

TaqMan® Copy Number Assays include pre-designed collections for both human and mouse genomes. The human collection includes more than 1.6 million pre-designed assays with genome-wide coverage targeting: known genes (exons, introns, and junctions), Copy Number Variant (CNV) sequences within the Database of Genomic Variants (DGV), and extragenic/non-gene regions. The mouse collection includes more than 180,000 assays targeting gene exons. Pre-designed assays to common vector marker and reporter genes are also available for transgenic studies. Pre-designed assays have been checked for genome quality as well as reference assay compatibility.



Figure 1. Determining copy number using CopyCaller™ Software.



TaqMan® Copy Number Assays, Custom Plus TaqMan® Copy Number Assays, and TaqMan® Custom Copy Number Assays.

Custom Plus TaqMan® Copy Number Assays are the optimal solution for studying variation in human and mouse genomic regions of interest for which a pre-designed assay is not available. A target range is defined on the Genome Map, then pre-masked targets are created and submitted to the assay design pipeline. Benefits include genome quality checks and human/mouse reference assay compatibility checks.

Custom TaqMan® Copy Number Assays are an option for additional targets of interest. The GeneAssist™ Copy Number Assay Tool enables you to submit your own pre-masked custom target sequences for assay design or primer/probe pair sequences for assay formulation. Custom assay designs do not go through genome quality checks, but can be compared with the human/mouse reference assays for compatibility in duplex reactions.

3. Quality

Whether testing a large or small number of targets against a few or hundreds of samples, having confidence in the data you produce is paramount. The Copy Number Assay Design Pipeline has been optimized to create high-performing copy number assays based on R&D wet lab validation of assays. Additionally, Life Technologies has developed the TaqMan® Assays QPCR Guarantee to provide you with confidence in the data you generate with a pre-designed TaqMan® Assay. Life Technologies guarantees that TaqMan® Copy Number Assays will perform to your satisfaction. If you are not satisfied with the performance of a TaqMan® Assay, we'll replace it at no cost, or credit the customer account for the purchase price of the assay.* Learn more at www.appliedbiosystems.com/taqmanguarantee.

4. Convenience

Choose only the specific assays you need for your project from the large variety of pre-designed assays, select a genomic target region of interest for Custom Plus Assay design, or submit your target sequence for Custom Assay design. The custom workflow also allows for the submission of primer/probe sequences for assay formulation.

TaqMan® Copy Number Assays are available in made-to-order small, medium, and large sizes.

TaqMan® Copy Number Assays (excluding Reference Assays) are shipped in a pre-formulated liquid format at ambient temperature and arrive ready to use (no need to resuspend different components). There is no optimization required; just use universal cycling conditions for every assay.

Ordering Method

TaqMan® Copy Number Assays and other related products can be ordered online at www.appliedbiosystems.com.

For help with online ordering, please refer to the Assay Search and GeneAssist™ Copy Number Assay Workflow Builder help text for Pre-designed, Custom Plus, and Custom Assays.

All TaqMan® Copy Number Assays are provided with a data CD containing:

- Protocol for the TaqMan® Copy Number Assays (PDF)
- Quick Reference Card for the TaqMan® Copy Number Assays (PDF)
- Protocol for CopyCaller™ Software (PDF)
- Quick Reference Card for CopyCaller™ Software (PDF)
- Product Insert for the TaqMan® Copy Number Assays (PDF)
- Assay Information File (AIF)
- Data Sheet (PDF)
- Safety Data Sheet(s) (PDF)
- Understanding Your Shipment Document (PDF)

Applied Biosystems® Instrument Compatibility

ViiA™ 7 Real-Time PCR System
 7900HT Fast Real-Time PCR System
 7500 Fast Real-Time PCR System
 7300 Real-Time PCR System
 StepOnePlus™ Real-Time PCR System

References and Additional Reading

Product Bulletin: TaqMan® Copy Number Assays (CO15338)
 Application Note: Design Pipeline and Validation of TaqMan® Copy Number Assays (135AP03-01)
 Protocol: TaqMan® Copy Number Assays (4397425)
 Quick Reference Card: TaqMan® Copy Number Assays (4397424)
 Quick Reference Card: CopyCaller™ Software (4400043)
 User Guide: CopyCaller™ Software (4400042)

Additional Product Information

Online Resources Page: www.appliedbiosystems.com/cnv

Copy Number Variation Web Resources

Database of Genomic Variants: <http://projects.tcag.ca/variation/>
 Sanger Institute Copy Number Variation Project: <http://www.sanger.ac.uk>
 UCSC Genome Bioinformatics Site: <http://genome.ucsc.edu/>
 Ensembl: <http://www.ensembl.org/>
 dbVar: <http://www.ncbi.nlm.nih.gov/dbvar/>

GeneAssist™ Copy Number Assay Workflow Builder

The GeneAssist™ Copy Number Assay Workflow Builder enables users to search for pre-designed human, mouse, and marker/reporters assays, to select a genomic human or mouse target region for Custom Plus Assay design, or to submit additional targets of interest for standard custom design/synthesis.

CopyCaller™ Software

CopyCaller™ Software was developed specifically for TaqMan® Copy Number Assay data analysis to enable users to obtain results quickly, robustly, and definitively (Figure 1). This free, easy-to use software utilizes a graphical interface and quickly calculates the raw and possible copy numbers for a set of samples. It also gives a confidence value for each copy number call and has outlier removal functionality. CopyCaller™ Software can perform relative quantitation analysis using a known calibrator sample or without a calibrator sample.




Figure 2. The complete TaqMan® Copy Number Assay workflow. TaqMan® Copy Number Assays have a simple workflow, and on average, setup to primary analysis takes only 3 to 4 hours (including a ~2 hour PCR run). Each copy number quantitation reaction contains four components: a TaqMan® Copy Number Assay, a TaqMan® Copy Number Reference Assay, TaqMan® Master Mix, and purified genomic DNA sample, run in four replicate wells.

Ordering Information

TaqMan® Copy Number Assays		Number of Reactions		Part Number		
Assay Scale	Concentration	384-well, 10 µL	96-well, 20 µL	Pre-designed Assays	Custom Plus Assays	Custom Assays
Small	20X	720	360	4400291	4442487	4400294
Medium	20X	1,500	750	4400292	4442520	4400295
Large	60X	5,800	2,900	4400293	4442488	4400296

TaqMan® Copy Number Reference Assays		Number of Reactions		Part Number
Description	Concentration	384-well, 10 µL	96-well, 20 µL	
Human Assays				
RNase P, 750 Rxns	20X (1 tube)	1,500	750	4403326
RNase P, 3,000 Rxns	20X (4 tubes)	6,000	3,000	4403328
TERT, 750 Rxns	20X (1 tube)	1,500	750	4403316
TERT, 3,000 Rxns	20X (4 tubes)	6,000	3,000	4403315

TaqMan® Copy Number Reference Assays		Number of Reactions		Part Number
Description	Concentration	384-well, 10 µL	96-well, 20 µL	
Mouse Assays				
Tfrc, 750 Rxns	20X (1 tube)	1,500	750	4458366
Tfrc, 3,000 Rxns	20X (4 tubes)	6,000	3,000	4458367
Tert, 750 Rxns	20X (1 tube)	1,500	750	4458368
Tert, 3,000 Rxns	20X (4 tubes)	6,000	3,000	4458369

*Subject to certain restrictions and terms and conditions. Visit www.appliedbiosystems.com/taqmanguarantee for details.

For Research Use Only. Not for use in human or animal therapeutic or diagnostic procedures.

NOTICE TO PURCHASER: LIMITED LICENSE

Practice of the patented 5' Nuclease Process requires a license from Applied Biosystems. The purchase of TaqMan® Copy Number Assays includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research when used with the separate purchase of an Authorized 5' Nuclease Core Kit. No other patent rights are conveyed expressly, by implication, or by estoppel. For further information on purchasing licenses contact the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94040, USA.

© Copyright 2010 Life Technologies Corporation. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan® is a registered trademark of Roche Molecular Systems, Inc. CO15424 1210

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in North America 800 955 6288

www.lifetechnologies.com

